

BIOCREDIT COVID-19 Ag

SARS-CoV 2 Antigen test

- » The first commercialization Black Gold Particle technology In the world (high sensitivity and specificity)
- » User friendly Dual Color System (Control line: Red, Test line: Black)
- » Rapid Test Time (5-8 Min)
- » Simple detection of SARS-CoV2 antigen

Core Technology Black Gold Particle



Differentiated Black Gold
Conjugate Technology



Dual Color System

8 minutes is ok!



BIOZOL
FIT FOR SCIENCE

BIOCREDIT COVID-19 Ag rapid test kit

Intended Use	Detection of SARS-CoV2 antigen 20
Package	Tests / kit
Storage	1~40°C
Specimen Type	Nasopharynx / Nasopharyngeal
Shelf life	12 months from manufacture
Time to result	date 5-8 minute

COVID-19 Rapid Diagnostic Test

PROCEDURE

Insert the swab specimen and swirl the swab 5-10 times.

Remove the swab while gently squeezing the head of the swab.

Close the assay diluent tube with a filter cap securely.

Invert the assay diluent tube and gently squeeze it to draw 3-4 drops (90-150µl) into a sample well on the device.

Read the result within 5-8 minutes.



INTERPRETATION

NEGATIVE	POSITIVE	INVALID
One red line "C" within the result window.	Two bands, black To test line and red "C" control line within the result window.	No "C" line within the result window. It is recommended that the specimen be retested.

BIOZOL
FIT FOR SCIENCE

BIOZOL ist
ISO 9001:2015
zertifiziert

BIOZOL Diagnostica Vertrieb GmbH

Leipziger Straße 4, 85386 Eching, Germany
T.: +49 (0)89 3799666-6 | info@biozol.de

www.biozol.de



RapiGEN RapiGEN INC.